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Instructions for use

SKY temp / Temporary restorations

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

The SKY temp abutments are prefabricated abutments made of plastic to be used for the SKY classic and blueSKY implant systems. They consist of a plastic abutment and a locking screw and are available in shoulder (neck) sizes S, M and L.

2. Indication / Application areas

The SKY temp abutments are used as base abutments for temporary restorations made of resin/composite but can also be used as gingiva formers.

3. Contraindication

The SKY implant system is not suitable for standard contraindications in implant dentistry such as patients with cardiovascular problems, bruxism, osteoporosis, heavy smokers, alcoholics, diabetes and metabolic diseases and insufficient primary stability of the implants.

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

SKY temp abutments should not be used for definitive prosthetic restorations, long-term temporaries (with a period of wearing of more than 6 months) and bridge restorations with more than 2 bridge units between the implants.

4. Storage and durability

The product is supplied in the non-sterile condition and dry and dust-free storage in the original packaging must be ensured.

5. Processing / Use

Generally, we recommend thorough surgical planning and selection of the prosthetic components depending on the situation and type of restoration.

The retention screws used in the laboratory are not suitable for the clinical use, since they are subject to stress in the laboratory; as a result long-term and reliable fit can no longer be guaranteed. bredent medical offers special laboratory screws to be used in the dental laboratory.

All components must be disinfected prior to the insertion (see item 8).

The SKY temp abutments are selected according to the desired gingiva profile and the width of the gap. The shape, emergence profile and height of the abutments can be customized. They can be adapted to the oral situation using the standard treatment methods and are screwed onto the implant and tightened with a maximum torque of 18 Ncm. Temporary crowns or bridges are fastened on the SKY temp abutments using temporary cement.

If the SKY temp abutments are used as gingiva formers, they need to be selected in accordance with the required gingival profile and the width of the gap. The shape, emergence profile and height of the abutments can be customized and they need to be screwed onto the implant and tightened with a maximum torque of 18 Ncm.

6. Technical data

REF	Product	Material
	SKY temp	
SKYTEMPS	SKY temp S incl. screw 2,2 mm	POM (Polyoxymethylene)
SKYTEMPM	SKY temp M incl. screw 2,2 mm	POM (Polyoxymethylene)
SKYTEMPL	SKY temp L incl. screw 2,2 mm	POM (Polyoxymethylene)
	Screw 2,2 mm	Titanium grade 4
	SKY Prosthetics accessories	
SKYLPS22	SKY laboratory screw 2.2 grey	Titanium grade 4, grey

7. Symbols

	Lot number
	Observe instructions for use
	Single use
	Non-sterile
	Dry storage

8. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff.

Only original tools and components may be used for processing.

The lot numbers of all components should be recorded to ensure traceability and to handle complaints.

Prior to the integration into the patient's mouth, all components must be disinfected. Only standard disinfectants for dental instruments may be used, e.g. Dürr ID 212 forte. The manufacturer's instructions on the use of the disinfectant must be observed.

After disinfecting, disinfectant residues must be thoroughly removed from the components (see the instructions of the manufacturer of the disinfectant; for example under running water - for ID212 forte).

All products labelled  are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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